

# SABIÖ 140

Automated ESR Analyzer  
For the direct determination from blood/EDTA Samples



# SABIO 140

## Automated ESR Analyzer

SABIO140 Automated ESR Analyzer is a bench top system which has complete integration in the hematology laboratories by the direct determination of ESR from blood/EDTA samples with its smart and innovative technology combining “ccd cam” and internationally accepted Westergren method.

SABIO140 is completely cost effective, provides patient comfort and environment friendly since it is using same EDTA tubes coming from CBC tests.

40 sample loading positions, suitability to continuous loading for 23 minutes when you start running the instrument, being maintenance free, being environment-friendly, LCD touch screen and data transmission to LIS are just some of the advanced features of the analyzer which provides laboratory comfort and productivity.

- **Use of the same CBC tubes. No dedicated citrate tube needed**
- **Westergren method by CCD camera**
- **Use of smart cards**
- **No reagents required, environmentally friendly**
- **40 sample loading positions, up to 120 tests/ hour**
- **Continuous loading of samples until the first result**
- **First result available after 24 minutes and each following result within 5 sec.**
- **10 inch LCD touch screen**
- **User friendly software**
- **Built-in thermal Printer (optional external printer)**
- **LIS connection**
- **Optional Barcode reader**

# RELIABILITY | PRODUCTIVITY | SATISFACTION




## USE OF CBC TUBES

No dedicated ESR citrate tube needed.  
No reagents required.



## EFFICIENT & USER FRIENDLY

Easy to use. 40 sample positions  
and continuous loading



## FEEL THE EXCELLENCE...

# SABIO 140

Automated ESR Analyzer

## Technical Specifications

Principles	: Westergren method by CCD Camera
Sample Loading Positions	: 40, Direct loading of Blood/EDTA Samples
Throughput	: Up to 120 samples/hour
Time to First Result	: First result in 24 minutes
Time to Result	: 5 seconds
Minimum Sample Volume	: 2 ml
Sample Test Volume	: 2 ml-3 ml
Barcode Reader	: External (optional)
Internal Printer	: Yes
LIS Communication	: Yes
Display	: Liquid Touch Screen , 800x600 pixel
Operating Environment	: Temperature: 15°C-35°C, Humidity: ≤80%
Power Requirement	: A.C. 110-230V ≤300 VA 50/60Hz
Dimension	: 690mm(L) x 425mm(W) x 315mm(H)
Weight	: 28KG